

ADVAPIX QUAD Version 1.0 - Datasheet

WITH CENTRAL OPENING











Illustration of single particle sensitivity of Timepix device. The tracks of different particles of radiation background (mostly muons and few protons) were recorded in 5 minutes on board of airplane. No noise (clean zero) is seen in dark regions. Inset shows the time profile along one muon track. The **ADVAPIX QUAD** camera contains four single Timepix devices with fast parallel readout with a 4 mm² central opening to allow the beam pass through the device uninterrupted. The AdvaPIX QUAD allows a sustainable frame rate of up to 1300 frames per second. Each Timepix device has an independent USB 3.0 communication channel. That assures fast read-out of the complete system. All modules in the system are operated synchronously or triggered independently. The sensor type and thickness can be customised.

The fast modules with Si or CdTe pixel detectors Timepix can be also used in different configurations such as stack of many layers or tiling to cover larger area or combination of both. Each module contains single Timepix device with fast parallel readout up to per second independent of frame occupancy. A separate USB 3.0 communication channel for each device assures fast read-out of the whole modular system. All modules in the system can be operated synchronously or triggered independently. The sensor type and thickness is of customer's choice.

The **ADVAPIX QUAD** is specially designed for combined Wide Angle X-ray Scattering (WAXS) and Small Angle X-ray Scattering (SAXS). The open space in the center of the camera eliminates completely the need to use a beam stop in front of the camera.

Main Features:

- Readout chip type..... Timepix
- Pixel size55 x 55 μm
- Sensor resolution......256 x 256 pixels
- Dynamic range in one frame......11 810¹
- Sensor material......100,300,500 μm Si or 1000 μm CdTe
- Power..... External or via second USB 3.0
- Interface.....USB 3.0 (Super-Speed)
 - Maximum readout speed......1300 fps
- Weight 1500 g

¹ A final picture can be created as sum of individual images. The dynamic range is limited only by the measurement time.





Device Parameters

Operating Conditions

Symbol	Parameter	Min	Тур	Max	Units	Comment
ТА	Temperature Range	0	50	70	°C	
Φ	Humidity			60	%	Not condensing

Electrical Specification

 T_A = 25°C, USB voltage V_{CC} = 4.8V

Symbol	Parameter	Min	Тур	Max	Units	Comment
Vcc	Supply Voltage	28	24	20	V	
Icc	Supply Current					
Icc1	Chip active			1000	mA	
P 1	Power Dissipation			24	W	
BNC Inputs						
VINL	Voltage Low	0		0.7	V	
VINH	Voltage High	1.7		5.0	V	
BNC Outputs						
VOUTL	Output Voltage Low	0		0	V	
VOUTH	Output Voltage High		4.5	5.0	V	
Bias Voltage Source for Sensor Diode						
VBIAS	Bias Voltage	0		±450	V	Polarity is sensor dependent

Performance Characteristics for Timepix

Symbol	Parameter	Min	Тур	Max	Units	Comment
f	Frame Rate			1300	fps	with USB 3.0 cable
T _{READ}	Frame Readout Time ²		588		μs	with USB 3.0 cable
dT	Time resolution	20	100		ns	

² During Readout time (or Dead time), no charge is collected from the sensor.



Sensor parameters

T_A = 25°C

Symbol	Parameter	Si			CdTe	Units	Comment	
	Thickness	100	300	500	675	1000	μm	
σ	Energy resolution of energy discrimination threshold (σ @ 23 keV)	0.5		1.1	keV			
σ	Energy resolution of energy discrimination threshold (σ @ 6o keV)	0.6			1.5	keV		
σ	Energy resolution in full spectral mode (σ @ 23 keV)	0.7		1.3	keV			
σ	Energy resolution in full spectral mode (σ @ 60 keV)	1.0		1.5	keV			
	Typical detectable energy range for X-rays ³	5 to 60		5 to 600	keV	See chart below		
	Pixel size		5	5		55	μm	

X-ray attenuation in Silicon and CdTe



•0.3 mm of Silicon •0.5 mm of Silicon •1 mm of Silicon •1 mm of CdTe

³The maximum detectable energy is limited only by sensitivity of the selected sensor for the given radiation type. The maximum measured energy can reach several MeV in case of heavy charged particles





+24VDC connector

Main power supply (via standard 2.1 mm barrel connector).

USB connector

USB 3.0 type micro B, Standard USB 3.0 Super-Speed. All USBs must be connected to the computer.

BNC synchronization connectors

All the BNC input signals (if not used) are tight to its neutral state, where the offered feature is disabled.

- Ready In input signal (gate-in signal that disables the measurement when at logical zero).
- Ready Out equals one if the **ready in** signal is in logical one and at the same time the Device has finished current cycle (is ready for a new cycle). This signal may be used as **Busy out signal** commonly used in other systems.
- Trigger In initiates shutter according to **TriggerStg** value. Use only when Master not connected.
- Trigger Out mirroring shutter of the Device.

ΛΟΥΑΓΑΜ

Imaging the Unseen

Synchronization waveform example

Master signals	When all ready, Master initiates shutters	Shutter Time (may differ for master and slave according to settings)	Master reads data and when finished, asks slave if it is finished too	When ready, slave acknowledges his readiness	When all ready, Master initiates shutters	
Ready Out		20ns				
Trigger Out (Master Shutter)						
Ready In		200ns		/		
Slave Shutter		160ns	_			
Timing example		10ms	20ms		300ns	
	cycle N-1	Cycle N				

Note1: the decisive leading edges are indicated in orange color. Note2: WidePIX 2x2 is Master, cooperating device is Slave

Depending on the TriggerStg value, the behavior of shutter signal is following (see figure below for graphical visualization of described behavior).

TriggerStg value	Detector performs measurement (activated shutter), when:
o (log o)	Trigger-In signal is logical zero => the camera starts measurement and keeps measuring frames as long as the trigger-in signal is in logical zero.
1 (log 1)	Trigger-In signal is logical one => the camera starts measurement and keeps measuring frames as long as the trigger-in signal is in logical one.
2 (rising edge)	Trigger-In signal is changing from zero to one => the camera measures only one frame.
3 (falling edge)	Trigger-In signal is changing from one to zero => the camera measures only one frame.









Mechanical Dimensions

Without protecting cover



All dimensions are in mm.

Extreme care must be taken when removing protecting cover and handling the **ADVAPIX QUAD** without the protecting cover. Warranty does not apply to mechanical damage of the sensor and wirebonds.





With protecting cover



All dimensions are in mm.







Configuration Examples

WAXS/SAXS (Wide/Small angle X-ray scattering)

Combination of MiniPIX and AdvaPIX Quad - compact WAXS and SAXS system:







Configuration Examples

Single detector

(Single layer of 256 x 256 pixels, speed of 1700 frames per second)



Stacked detector

(Arbitrary number of layers of 256 x 256 pixels, speed of 1700 frames per second)



Quad detector or Quad detector with central hole

(Single layer of 512 x 512 pixels, speed of 1300 frames per second with optional hole in the middle, up to 2mm square)



Row detector (Single layer 256 x 1024 pixels, speed of 1700 frames per second)





Warning

Do not touch sensor surface!

Instructions for safe use. To avoid malfunction or damage to your **ADVAPIX QUAD** please observe the following:

• Do not expose to water, moisture.

• Do not disassemble. Wire-bonding connection may be irreversibly damaged.

• Do not insert any object into the sensor window.

Copyright





DC@WH97`;aV<`````Dc`mhYWD`Uhn`%!`+`````8`!+*'''+`KU`XVfcbb`````;9FA5BM` HY`.Ž(-`f#&(''Ł*\$('%+'\$````:UI.`Ž(-`f#&(''Ł*`--`((```9!AU]`.`cY_4 dc`mhYWXY````kkk"dc`mhYWXY

