



## Software Development Kit



SDK (Software Development Kit) provides the user with the detailed Function Calls to our .DLL files. This package is designed for customers who wish to create their own custom software interface allowing complete control over your spectrometer or system.

Every spectrometer / system we sell can be run using the SDK, including RS232 and USB units, from non-cooled spectrometers to complete Raman systems with Laser Power Control.



### SDK Program Includes:

The SDK Package comes complete with simple programming examples done with our various spectrometer models to get you started.

- C#
- C++ Builder 6
- Visual C++ 6.0
- Visual Basic 6.0
- VBA
- Labview
- VB.NET

### SDK Supports:

- Raman Systems
- Reflectance Systems
- Transmittance Systems
- USB Cleanlaze® Laser
- RS232 Spectrometers
- USB Spectrometers

### Program Examples:

Specification	Interface
C# SDK	USB Spectrometers
VB.Net SDK	USB Spectrometers
C++ Builder 6 SDK	USB Spectrometers
Visual C++ 6.0 SDK	USB Spectrometers
Visual Basic 6.0 SDK	USB Spectrometers / Lasers
VBA SDK	USB Spectrometers
Labview 8.2 SDK	USB Spectrometers
RS232 Interface SDK	RS232 Spectrometers

