

HSV-100 High Speed Vibrometer

The HSV-100 High Speed Vibrometer measures vibration displacements and velocities up to 40 m/s, simultaneously, with an arbitrary number of measurement channels.

Vibration of the sample can be measured reliably even when space is restricted, thanks to the very compact optical sensor head coupled to a separate laser unit. One benefit of this non-contact method is that structural motion is measured more accurately.

State of the art digital decoding in the controller achieves high resolution signal outputs with low noise for clear and precise results.

Several dedicated HSV-100 models are available for different measurement tasks: The single-channel version is the solution for common high speed vibration measurement. The dual-channel configuration allows differential vibration measurement, which automatically compensates for any background vibration of the test rig. An arbitrary number of extra channels can be added for true multi-channel differential vibration analysis.



Benefits

- High resolution real-time measurement
- Intuitive “point and shoot” operation
- Two- and multi-channel measurement with correct phase relationships
- IP64 rated sensor head housing for extreme industrial applications
- Frequency bandwidth 50 kHz, 250 kHz optional
- Eye-safe class 2 laser

HSV-100 High Speed Vibrometer

Optical Vibration Measurement
at High Speeds

Datasheet



Technical Data



Controller Specifications

Controller model	HSV-E-100-01	HSV-E-100-02	HSV-E-100-01MC
No. of measurement channels	1: Channel A	2: Channel A, Channel B	1: Channel A ¹
Signal outputs ²	Velocity: A Displacement: A	Velocity: A, B, A-B Displacement: A, B, A-B	Velocity: A, B, A-B Displacement: A, B, A-B
Analog velocity output	BNC, +/-8 V	BNC, +/-8 V (+/-12 V for A-B)	BNC, +/-8 V (+/-12 V for A-B)
Analog displacement output	BNC, +/-10 V	BNC, +/-10 V	BNC, +/-10 V
Frequency bandwidth	0 Hz ... 50 kHz; 0 Hz ... 50 kHz / 250 kHz switchable (optional)		
Max. velocity ³	40 m/s	40 m/s	40 m/s
Velocity ranges	8 ranges: 25 mm/s/V ... 5 m/s/V		
Displacement ranges	11 ranges: 0.25 mm/V ... 256 mm/V		
Tracking filter	4 settings (slow, medium, fast, off)		
Signal level	Bargraph at front panel; Output as voltage signal (BNC, 0 ... 5 V)		
PC interface	USB 1.1, remote control of controller settings		

Optical Specifications

Sensor head model	HSV-I-100-FF	HSV-I-100-V25	HSV-I-100-V40
Laser	Helium Neon (HeNe), 633 nm, visible red laser beam		
Laser safety class	Class 2 (eye-safe), <1 mW		
Working distance	fixed: 520 mm	variable: 110 mm ... approx. 10 m	variable: 200 mm ... approx. 10 m

General Specifications

System component	HSV-E-100 Controller	HSV-I-100 Sensor Head: Laser Unit	HSV-I-100 Sensor Head: Sensor
Dimensions [L x W x H]	450 mm x 360 mm x 150 mm (19", 84HP/3U)	340 mm x 130 mm x 115 mm (w/o handle)	174 mm x 48 mm x 39 mm (HSV-I-100-FF Sensor)
Weight	9.5 kg	6 kg	0.5 kg
Protection class	IP20	IP64	IP64
Cable length	Electrical main cable between Controller and Sensor Head: 10 m (detachable); Fiber cable between Laser Unit and Sensor: 3 m (permanent)		
Operating temperature	+5 °C ... +40 °C		
Storage temperature	-10 °C ... +65 °C		
Relative humidity	max. 80%, non-condensing		
Mains voltage	100 ... 240 VAC +/-10%, 50/60 Hz		
Power consumption	max. 70 VA		

- Additional measurement channel for HSV-E-100-02 Dual Channel Controller; the HSV-E-100-02 can be extended by an arbitrary number of additional channels.
- All signal outputs are available simultaneously.
- For the differential signal A-B velocities up to 60 m/s are measurable due to the extended voltage range of the signal output.



Options and Accessories:

HSV-AK-800 Beam Steering Unit

Optical beam deflection system for easy alignment of the laser beam especially for valve train measurement applications; compatible with the HSV-I-100-FF Sensor Head.

Polytec GmbH (Germany)
Polytec-Platz 1-7
76337 Waldbronn
Tel. +49 7243 604-0
info@polytec.de

Polytec GmbH (Germany) Vertriebs- und Beratungsbüro
Schwarzschildstraße 1
12489 Berlin
Tel. +49 30 6392-5140

Polytec, Inc. (USA)
North American Headquarters
16400 Bake Parkway
Suites 150 & 200
Irvine, CA 92618
Tel. +1 949 943-3033
info@polytec.com

Central Office
1046 Baker Road
Dexter, MI 48130
Tel. +1 734 253-9428

East Coast Office
25 South Street, Suite A
Hopkinton, MA 01748
Tel. +1 508 417-1040

Polytec Ltd. (Great Britain)
Lambda House
Batford Mill
Harpenden, Herts AL5 5BZ
Tel. +44 1582 711670
info@polytec-ltd.co.uk

Polytec France S.A.S.
Bâtiment Orion – 1er étage
39, rue Louveau
92320 Châtillon
Tel. +33 1 496569-00
info@polytec.fr

Polytec Japan
Arena Tower, 13th floor
3-1-9, Shinyokohama
Kohoku-ku, Yokohama-shi
Kanagawa 222-0033
Tel. +81 45 478-6980
info@polytec.co.jp

Polytec South-East Asia Pte Ltd
Blk 4010 Ang Mo Kio Ave 10
#06-06 TechPlace 1
Singapore 569626
Tel. +65 64510886
info@polytec-sea.com

Polytec China Ltd.
Room 1026, Hanwei Plaza
No. 7 Guanghua Road
Chaoyang District
100004 Beijing
Tel. +86 10 65682591
info-cn@polytec.com